


# Communicating with External Hardware using PsychoPy

 by Kimberley Dundas



# What kind of external hardware?

- PsychoPy is currently able to communicate with:
- EEG recording hardware
- fMRI/MRI scanners
- Eye-trackers
- Microcontrollers (e.g. an [Arduino board](#))
- fNIRS – there's a really clear [video tutorial](#) from NIRx

# EEG recording hardware - caveats

- Timing issues and synchronisation
  - Timing Mega-study (Bridges et al., 2020)
  - Dropped frames
  - Millisecond precision
  - Use Builder!
  - Test, test, test:
    - Monitor refresh rates – Poth et al. (2018)
    - ‘Drawing’ from top to bottom
  - Lag and variability

# How do I send triggers/markers?

## ▼ Parallel port component

8 bit communication

Sends values in parallel - faster

Not all computers have a parallel port...

## ▼ Serial port component

Sends values serially – slower?

Adapters are available for USB, e.g.:

[LabHackers USB2TTL8](#)

[LabJack U3](#)

# Sending EEG triggers in PsychoPy

The biggest difficulty that's faced when trying to send triggers in PsychoPy is working out:


1. Which port address your EEG acquisition system is connected to
2. How your EEG acquisition system *wants* to receive the triggers you send

Once we know these things - sending triggers is relatively simple!



Be sure to bear the timing caveats in mind when designing your EEG experiments

# Serial Port Component



Insert a **Serial Out** component to the routine that you'd like to add triggers to

**Start** refers to when you want the trigger to be sent; this can be at a certain time or frame, or we can set the trigger to be sent when a certain condition is met (more on this later)


**Port** refers to the address of the serial port your device is connected to

**Start data** refers to the value that you want to actually send to your acquisition system - exactly what you'd like to send will depend on what your system wants to receive. This can take a variable just like most other fields so that you can send different triggers for different stimuli etc.


# Parallel Port Component

The parallel port component is essentially set up in the same way as the serial port component, but the fields are spread across separate tabs:


**Start** is set in the Basic tab just as with the serial port component:




**Start data** can be found in the Data tab



**Port address** is found in the Hardware tab




- ☐ If you don't see your port address listed, you can add it by navigating to File > Preferences > Hardware > Parallel Ports > ... > New:



# Sending triggers to mark events

Exactly *when* your trigger will be sent to your acquisition device will depend on how you set the **Start** of the serial or parallel port component:



We might want to send our trigger on a certain frame or time, but typically we want to send a trigger to mark a given event. To do this we have to set the Start to **Condition** (rather than Duration)

To send on the **onset of something** (e.g., a visual stimulus) set that condition to be (where stimulus is the name of the component you want to yoke your trigger to):

```
stimulus.status == STARTED
```

To send on a **keypress**, keep the Start as Condition and set the condition to be (where key\_response is the name of your keyboard component):

```
key_response.keys
```

To send when a stimulus is **clicked by a mouse** (mouse here refers to the name of your mouse component, and stimulus\_mouse refers to the stimulus you want to be clicked):

```
mouse.isPressedIn(stimulus_mouse)
```


# Using MRI/fMRI with PsychoPy

- Keypress
- Serial/parallel communication
- Find out more (including code to copy and paste) [here](#).

② Non-slip timing – what is it and why do we care?




# Relative timing



Set a timer to **0** at the start of each **trial**





Present a stimulus  
Check on each frame whether the timer exceeds our intended duration  
If it hasn't, present the stimulus for this frame too



Stop on the frame where the timer **exceeds** our intended duration

Potentially leading to an overshoot, for example of 10ms

# Non-slip timing



Set a timer that is **not** reset at  
the start of each **trial**

Instead, any residual value from  
the previous trial is added to the  
intended duration of this trial,  
e.g. -10ms

Present a stimulus

Check on each frame whether  
the timer exceeds our intended  
duration


If it hasn't, present the stimulus  
for this frame too

Stop on the frame where the  
timer **exceeds** our intended  
duration


If this was 10ms too long, our  
residual is -10ms, which is then  
added to the next trial

# Adding eye-tracking to your experiment

# Mouse Gaze




# Calibration and Validation routines




- ⚠️** These components create “standalone routines” rather than components. Once selected you select “insert routine” and add the calibration routine where you want it on your flow.

# Start the eye tracker




You can start/stop the eye tracker in a single trial – or start the tracker before your trials and stop it afterwards.

# Stop the eye tracker




Here we start the tracker in the instructions and stop it in the "bye" routine.

# Using ROI components




We can set our trial to end when the participant looks at the ROI for a set period of time...

# Using ROI components




Since we know where the target is on every trial, we can update the region of interest location on every trial too...

# Feedback based on ROI data



# Configure your experiment to your eye tracker



# Exercise (10 mins)

- ⓘ Make the picture of the magnifying glass track/move with your eye position.